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Grinding must be performed in an area separate from the testing area. Use the Romer Mill - Model 2A, Bunn Grinder, or equivalent to grind the sample.

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
## Illustration of Positive and Negative Results

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## NON GMO TESTING QUICK COMB TESTING

**RED--ROUND UP READY** -- Retest allowed, if fails 2nd test  
corn can go as regular.

**GREEN--B.T.--** Retest allowed, if fails 2nd test corn can go  
as regular.

**BLUE--STARLINK--** NO RE TEST ALLOWED. Can not accept corn at all.

**ORANGE--YIELD GARD : ROOTWORM--** Can go as regular if fails.  
NO RETEST

**PURPLE--CRY1F--** NO RETEST CORN CAN BE DUMPED AS REGULAR. -Herculex-

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
**1: Test Portion** (Original testing service)

**2: File Portion** (Review testing service)

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
		Grain Inspection, Packers and Stockyards Administration	P.O. Box 640 Destrehan, Louisiana 70047 (504)-764-2326
SUBJECT: Confirmation of Telephone Report			
TO: Food and Drug Administration New Orleans District Office			
This will confirm our telephone report of ____ \ ____ \ ____ regarding the following Export lot of ____ (Product) <b>Temporary FDA Fax # (504) 219 - 8813</b>			
FDA Contact person: Ms. Patricia Schafer <b>Temporary FDA Phone # (504) 219 - 8818 ext 101</b>			
Elevator Name: _____		Carrier Identification: _____	
Commodity: _____		Destination: _____ Sampling Date: _____	
Elevator Contract #: _____		FGIS Load Order #: _____	
Original Inspection results: _____ PPB - REX result: _____ Test method: _____			
Appeal results: _____ PPB - Test method: _____			
Reconditioned results Whole Corn: _____ PPB; Reconditioned Screenings _____ PPB			
<b>Export Lot Disposition</b>			
Identification of Grain	Amount LBS	Lot Disposition Sealed in Bin No.	Remarks
Shipping Bin Nos. _____ Orig wt. _____	Transfer to secured H/B weight: _____	Bin# _____ Seal # _____ Bin# _____ Seal # _____ Bin Top Seal # _____	Checked empty by _____ Time: _____ Date: _____
<b>Reconditioning Whole Grain or Transfer</b>	Transfer to secured H/B weight: _____	Bin# _____ Seal # _____ Bin# _____ Seal # _____ Bin Top Seal # _____	Checked empty by _____ Time: _____ Date: _____
<b>Reconditioned Screenings or Transfer</b>	Transfer to secured H/B weight: _____	Bin# _____ Seal # _____ Bin# _____ Seal # _____ Tank# _____ Seal # _____	Checked empty by _____ Time: _____ Date: _____
Original results Faxed to FDA on ____ \ ____ \ ____ By _____ Elevator Ph # _____			
Reconditioning results Faxed on ____ \ ____ \ ____ By _____ Ph# _____			
Remarks: FDA File Sample pick-up by: _____ Date: _____			
Revised 05-30-08:		See reverse side for additional flow path documentation (original copy only).	

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AFLATOXIN HANDBOOK  
CHAPTER 1  
12-18-06

## 1.5 APPROVED TEST METHODS

FGIS has approved test kits for use at field testing locations. The Aflacard T20, AflaCup, AgraStrip, Agri-Screen, QuickTox, Reveal, and ROSA® P/N test kits are approved for qualitative analysis of corn. The Aflatest, Fluoroquant, Myco✓, RIDASCREEN Fast Aflatoxin SC, ROSA® Quantitative, Veratox, and Veratox-AST test kits provide quantitative analysis but can be used for qualitative results. High Performance Liquid Chromatography (HPLC) testing is reserved for quantitative testing at the Technical Services Division (TSD) only.

FGIS APPROVED TEST METHODS			
Method and Test Kit	Approved for		Test Kit Range
	Qualitative	Quantitative	
Aflacard T20 – (R-Biopharm Rhone)	X		20 ppb
AgraStrip – (Romer)	X		10 or 20 ppb
AgriScreen – (Neogen)	X		20 ppb
QuickTox – (Envirologix)	X		20 ppb
Reveal for Aflatoxin – (Neogen)	X		20 ppb
ROSA® Aflatoxin P/N – (Charm Sciences, Inc.)	X		10 or 20 ppb
Aflatest – (Vicam)	X	X	5 - 100 ppb
Fluoroquant – (Romer)	X	X	5 - 100 ppb
Myco✓ – (Strategic Diagnostics Inc.)	X	X	5 - 80 ppb
RIDASCREEN Fast Aflatoxin SC – (R-Biopharm)	X	X	5 - 100 ppb
ROSA® Aflatoxin Quantitative – (Charm Sciences, Inc.)	X	X	5 - 50 ppb
Veratox – Neogen	X	X	5 - 50 ppb
Veratox AST – (Neogen)	X	X	5 - 300 ppb

**NOTE:** The test ranges are for performing an individual analysis with an undiluted sample extract. To obtain accurate results above the test kit range a supplemental analysis must be performed.

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Information For...

Browse By Subject

- ▶ Grain, Rice & Pulses
- ▶ Livestock, Meat & Poultry
- ▶ Education & Outreach
- ▶ Laws & Regulations
- ▶ Markets & Trade

You are here: [GIPSA Home](#) /

## News & Events

### GIPSA APPROVES R-BIOPHARM RIDASCREEN® FAST SC AFLATOXIN SC TEST KIT TO DETECT AFLATOXIN IN RICE BRAN

WASHINGTON, Aug. 03, 2007 - The U.S. Department of Agriculture's Grain Inspection, Packers and Stockyards Administration (GIPSA) today approved the R-Biopharm, Inc., Ridascreeen® FAST Aflatoxin SC quantitative test kit, Product # R6001, for the official testing of rice bran for total aflatoxins in the national grain inspection system.

"Approval to test additional commodities using this kit expands the choice of aflatoxin test kits available to the national grain inspection system," said GIPSA Administrator James E. Link.

The kit was previously approved for use on corn, cracked corn, corn flour, corn germ meal, corn gluten meal, corn meal, corn soy blend, distillers dried grains with solubles, popcorn, sorghum, soybeans, and wheat.

Aflatoxins are a toxin produced by a mold and can be present in corn, sorghum, and other grains.

Technical information: Lynn Polston, GIPSA Technical Services Division, (816)891-0444 or email [Lynn.A.Polston@usda.gov](mailto:Lynn.A.Polston@usda.gov).

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